



The spring has a rubber gromet that attaches to the handle. The spring and rubber gromet are not intended to stop the handle when released. The user should always use their hand to slow the release of the handle.

# Spring and Grommet Reinstallation



If the spring becomes dislodged using a small screwdriver it can easily be put back into place. Ensure the small end of grommet is exposed from the spring (Picture A). Then guide it into the handle (Picture B).

It's in our DNA<sup>®</sup>



# Spring and Gromet Reinstallation



Press the gromet and spring into the handle pocket. Once in place close the handle to ensure complete seating of gromet.



It's in our DNA<sup>®</sup>



# Completed Installation



The spring is now secured back to the handle.

